RESEARCH NOTE



Cryopreservation and validation of differentiated PDGFRα-positive cells for long term usage in experimentation

Ankita Srivastava¹, Reeja Rajan¹, Sallam Hasan Abdallah^{2,3}, Amir Ali Khan^{4,5} and Bashair M. Mussa^{6*}

Abstract

Objective Differentiation of immortalized Mesenchymal Stromal Cells (iMSCs) into PDGFRa-positive cells under controlled growth conditions has several vital implications in functional studies concerned with the pathogenesis of Diabetic Gastroparesis (DGP). A study published previously by our research group demonstrated the importance of these cells as a novel, in-vitro model for investigating the functional role of neuronal nitric oxide synthase. The currently available methods require fresh differentiation of PDGFRa-positive cells for each round of experimentation. This leads to longer delays, higher usage of reagents, and inconsistency in reproducibility of experiments frequently. We thus aimed to establish through validation that cryopreserving and maintaining the iMSC-derived PDGFRa-positive cells for functional investigations help us to overcome these challenges.

Results We demonstrated for the first time that the differentiated PDGFRa-positive cells from iMSCs can be cryopreserved and thawed to be used as per the experimental requirements with prolonged preservation of their characteristics. We assessed the viability of differentiated PDGFRa-positive cells pre- and post-freezing with the subsequent validation of their functional features using flow cytometry, qRT-PCR, and western blotting. We have been successful in demonstrating for the first time that the cryopreservation of previously differentiated PDGFRa-positive cells can be used as a feasible and cost-effective model for experimental reproducibility in functional studies of Diabetes Gastroparesis.

Keywords iMSCs, Stem cell differentiation, PDGFRa-positive cells, qRT-PCR, Flow cytometry, Western blotting

*Correspondence: Bashair M. Mussa bmussa@sharjah.ac.ae Full list of author information is available at the end of the article



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/A.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.



Introduction

A study published by our research group [1] has established that telomerase transformed immortalized Mesenchymal Stromal Cells (iMSCs) can be successfully differentiated into PDGFR α -positive cells using connective tissue growth factor (CTGF) and L-ascorbic acid (LAA) [1, 2]. The iMSC-derived PDGFR α -positive cells thus obtained have been successfully used as an in-vitro model to conduct functional studies associated with Diabetic Gastroparesis [1, 3]. We and others have used this protocol to differentiate iMSCs into PDGFR α -positive cells; however, one of the limitations of the protocol was the waiting time of 21 days for each round of differentiation. Therefore, we worked towards the idea of using cryopreservation as a method for maintaining these cells post day 21 of differentiation and the possible usage of these cells post thawing preferably after 3–4 weeks to render them useful in future experiments. This method would help investigators have ready to use cryopreserved differentiated PDGFR α -positive cells instead of having repeated rounds of 21-day differentiation to execute their experimentation. Figure 1 illustrates an overview of the differentiation, maintenance and cryopreservation process.

We have also described the details of the various steps used in validating the PDGFR α -positive cells during long-term storage in liquid nitrogen and their subsequent thawing and compared the results to a positive control (primary fibroblasts [F180]) and negative control (undifferentiated iMSCs). These steps involve analyzing cell survivability using MTT, detecting specific cell surface markers using flow cytometry, and analyzing gene and



Fig. 1 Cryopreservation, thawing and maintenance of differentiated PDGFRα-positive cells. iMSCs were treated with connective tissue growth factor (CTGF) and L-ascorbic acid (LAA) for a period of 21 days to allow their differentiation into PDGFRα-positive cells. Post day 21, the cells were cryopreserved for 3–4 weeks as a new approach, then thawed and sub-cultured. The survival of the cells was assessed using MTT assay and functional validation was conducted using flow cytometry for detection of cell surface markers, qRT-PCR for gene expression analysis and western blotting for protein expression analysis

protein expression using qRT-PCR and western blotting, respectively. The preservation of morphology of the differentiated cells were validated using Scanning electron microscopy.

Main text

Experimental design and methodology *Cell culture and fibroblastic differentiation*

(i) Reagents and consumables required

a. Negative control: telomerase transformed immortalized human bone marrow-derived mes-

enchymal stromal cells (iMSCs; RRID: CVCL_ B5PE); Applied Biological Materials, Richmond, Canada; Catalogue No. abm T0529.

- b. Positive control: primary human fibroblasts [F180]; University Hospital Sharjah, Sharjah, UAE.
- c. 6 well plates; Jet Biofil, Guangzhou, China.
- Minimum Essential Medium-alpha modification (MEM-α); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. 32571-028.
- e. Fetal Bovine Serum; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. F7524.
- f. L-Glutamine; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. G3126-500G.

- g. Penicillin–Streptomycin; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. P4333-100ML.
- h. Connective tissue growth factor (CTGF); Biovendor, Brno, Czech Republic; Ref Code: RD191035299R.
- i. L-Ascorbic acid (LAA); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. A92902.
- j. T-25 and T-75 cell culture flasks; Sigma Aldrich, St. Louis, MO, USA.
- k. Trypsin–EDTA solution; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. T4299.
- Dulbecco's Phosphate Buffered Saline (PBS); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D8537
- m. Dulbecco's Modified Eagle Media (DMEM); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D5671.
- n. Media recipes

1) Complete media: MEM-α and DMEM

To MEM- α (reagents #1.1d), 0.292 g/L of L-glutamine [reagents #3.1(i)f], 10% heat inactivated FBS [reagents #3.1(i)e] and 1% Penicillin–Streptomycin [reagents #3.1(i)g] were added. To DMEM [reagents #3.1(i)m], 10% heat inactivated FBS and 1% Penicillin–Streptomycin were added.

 Fibroblastic differentiation medium Complete MEM-α supplemented with 100 ng/L CTGF [reagents# 3.1(i)h] and 50 µg/ mL LAA [reagents #3.1(i)i].

(ii) Methodology

- 1) Culture telomerase transformed iMSCs [Reagents 3.1(i)a] in complete MEM- α (Media recipes n#1) in a humidified atmosphere of 5% CO2 at 37 °C in T-75 tissue culture flasks [Reagents #3.1(i)j].
- 2) Harvest the cells at 80% confluency by washing the adherent cells with 5 mL of sterile PBS thrice and trypsinize them with the addition of 500ul of 0.25% trypsin–EDTA solution [Reagents #3.1(i)k] and incubating for 2 min at 37 °C. Terminate the tripsinization process by the addition of 5 mL of complete DMEM [Media Recipes n#1].
- Plate the iMSCs in 6-well plates [Reagents #3.1(i)1b] at a density of 250×10³ cells/well.
- 4) Add fibroblastic differentiation medium when the cells reach 70–80% confluency (Media recipes n#2).
- 5) Change the medium every day for 21 days.

- 6) Use iMSCs [Reagents #3.1(i)a] cultured in MEM- α as negative control and primary human fibroblasts [F180] [Reagents #3.1(i)] cultured in DMEM as a positive control in the absence of CTGF and LAA.
- Subsequently, harvest the differentiated PDGFRα-positive cells for future studies on the 21st day for storage through cryopreservation.

Cryopreservation of the differentiated PDGFRa-positive cells

(i) Reagents and consumables required

- a. 1.2 mL cryovials; Thermo Fisher Scientific, Waltham, MA, USA; Catalogue No. 5000-0012.
- b. Trypsin–EDTA solution; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. T4299.
- c. Dulbecco's Phosphate Buffered Saline (PBS); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D8537.
- d. Complete MEM- α and DMEM media (Media recipes n #1).
- e. Dimethyl sulfoxide (DMSO); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D8418.

(ii) Methodology

- On day 21, wash the differentiated PDGFRαpositive cells twice with PBS [Reagents #3.2(i)c] to remove the old media.
- 2) Subsequently, add 500 μL of trypsin–EDTA [Reagents #3.2(i)b] and distribute over the entire surface area of the wells, ensuring complete coverage of the adherent cell layer and incubate at 37 °C for 2 min. Terminate the trypsinization process by the addition of complete DMEM [Media recipes n#1]
- Observe the cells under a microscope to confirm complete detachment, add fresh complete media (MEM-α, Media recipes n #1), disperse by gentle repeated pipetting to avoid clumping and achieve single-cell suspension. Similarly, for the controls, resuspend iMSCs in MEM-α and fibroblasts in DMEM (Media recipes n#1).
- 4) Collect the cell suspension in 15 mL falcon tubes and centrifuge at 1500 rpm (g-101) at 24 °C for 5 min.
- 5) CRITICAL STEP. Discard the supernatant and resuspend the pellets in 90% of fresh complete media (Media recipes n#1). To this, add 10% of cryoprotectant (DMSO; [Reagents #3.2(i)e] drop-wise to reduce cell death.

- Distribute the above cell suspension into cryovials with approximately 1×10⁶ cells per vial [Reagents #3.2(i)a].
- 7) Store the cryovials in a -80 °C freezer for 24 h up to 1 week and then transfer them to a liquid nitrogen tank until future use.

Thawing and sub-culturing of the differentiated PDGFRa-positive cells

(i) Methodology

- Thaw the frozen cell suspension in cryovial stored at - 80 °C/liquid nitrogen in a 37 °C water bath for less than a minute with gentle swirling. Transfer it into a 15 mL sterile tube and add 4–6 mL of complete media (Media Recipes n#1).
- 2) Centrifuge the above tube with cells at 1500 rpm for 5 min at 24 °C to get rid of DMSO.
- Discard the supernatant and resuspend the pellet in 4–6 mL of complete media. Transfer the cell suspension into a T-25 culture flask [Reagents #3.1(i) j]and incubate at 37 °C with 5% CO₂.
- 4) Over the following few days, observe the cells and change the media every two days until the cells are healthy and reach 60–80% confluency (up to 1 week).

Subculture the cells for one passage before using for further experimentation to achieve healthy differentiated cells. Change the medium (with CTGF and LAA) every 48 h.

Measuring cell viability using MTT

(i) Reagents and consumables required

- a) 3-[4,5-Dimethylthiazol-2-yl]-2,5 diphenyltetrazolium bromide (MTT); Sigma-Aldrich, St.Louis, MO, USA; Catalogue No. CAS 298-93-1.
- b) 96-well plates; Jet Biofil, Guangzhou, China.
- c) Trypsin–EDTA solution; Sigma Aldrich, St. Louis, MO, USA; Catalogue No. T4299.
- d) Dulbecco's Phosphate Buffered Saline (PBS); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D8537.
- e) Complete MEM- α and DMEM media (Media recipes n#1).
- f) Dimethyl sulfoxide (DMSO); Sigma Aldrich, St. Louis, MO, USA; Catalogue No. D8418.

(ii) Methodology

- 24 °C to get rid of DMSO.
 2) Seed approximately 1×10⁴ PDGFRα-positive cells in triplicates in 96-well plates [Reagents #3.4(i)b] and maintain in 200 µL complete MEM-α medium (Media recipes n#1) at 37 °C with 5% CO₂.
- 3) After 72 h, discard the old culture medium and add 20 μ L of MTT (5 mg/mL) [Reagents #3.4(i) a] dissolved in 100 μ L of PBS [Reagents #4.1d] to each well.
- 4) Discard the MTT+PBS solution carefully and add 100 μ L of DMSO; [Reagents #3.4(i)f] per well.
- 5) Record the absorbance at 570 nm on a microplate reader.

Results

Investigation of viability of cells post cryopreservation using MTT assay

The comparison of cell viability before and after cryopreservation of cells was done using the MTT assay and only a minute difference of 2.3% was observed in the recovered percentage of cells post cryopreservation, the difference being non-significant proves that the viability of cells is not much affected due to the cryopreservation process (Fig. 2).



Fig. 2 Comparing cell viability before and post cryopreservation. Using MTT cell viability assay, absorbance was recorded for assessing cell survival in group 1 (before freezing; 0.177 nm) and group 2 (after freezing; 0.154 nm). A minute difference of 2.3% was observed, which was not significant. Values represent average absorbance from three individual experiments

Confirming the presence of cell surface markers of differentiated PDGFRa-positive cells using flow cytometry analysis

The validation of the cryopreserved differentiated PDGFR α -positive cells was done by flowcytometry analysis. The expression of cell surface markers CD140 α /PDGFR α , CD44, and CD34 were measured and analysed in iMSCs, fibroblasts, and PDGFR α positive cells using flowcytometry. The details of the results can be obtained from published data [1] [https://doi.org/10.3390/ijms22073514].

Gene expression analysis for validation of the cryopreserved and differentiated PDGFRα-positive cells

The validation of the cryopreserved differentiated PDGFR α -positive cells was done by measuring the Gene expression of Extracellular Matrix (ECM) Proteins (COL I, DEC, ELA, HAS3, and TIMP1) and Stem Cell Differentiation (SCD) markers (ALP, AGG, CD44, and FSP-1) in fibroblasts, iMSCs, and PDGFR α -positive cells. The details of the results from qRT PCR can be obtained from previously published data [1] [https://doi.org/10.3390/ijms22073514] of our group.

Protein expression analysis of differentiated PDGFRa-positive cells using western blotting

The validation of the cryopreserved differentiated PDGFR α -positive cells was done by analysing the Protein expression of Stem Cell Differentiation (SCD) markers (ALP, AGG, CD44, and FSP-1) in fibroblasts, iMSCs, and differentiated PDGFR α -positive cells. The details of the results obtained from western blot analysis can be obtained from published data [1] [https://doi.org/10.3390/ijms22073514] (Additional file 1).

Morphological examination of the differentiated PDGFRa-positive cells post cryopreservation

The conservation in morphology of the cryopreserved differentiated PDGFR α -positive cells was studied using Scanning Electron Microscopy (SEM). The images of the morphological examination can be obtained from published data [1] [https://doi.org/10.3390/ijms220735 14].

Conclusion

The main advantages of developing the present protocol were reducing the time taken to differentiate PDGFR α -positive cells, making the cells readily available for functional studies, reducing the cost of reagents, and achieving an excellent reproducibility. MTT viability assay has revealed a non-significant difference in cell survival before and after freezing, indicating no adverse effects of long-term storage on the cells, which are comparable to the freshly differentiated cells. Validation experiments using flow cytometry have shown that the efficiency of differentiation is due to increased expression of the CD140a cell surface marker, also known as PDGFRa. Gene and protein expression analysis has also highlighted the intact functional capacity of these cells to be used for further experimentation. Comparing the outcomes with relevant positive and negative control cells cultured simultaneously further provides proof of maintained differentiation. However, future studies, including cell population analysis, are warranted to reveal important information related to the ratio of differentiated vs. non-differentiated cells in this homogenous population.

Limitations

One of the limitations of this study is obtaining high percentage of reproducibility between each round of cryopreservation and thawing. Therefore, each step of the differentiation process and thawing is critical. Secondly, although DMSO is considered as the best cryoprotectant for cell culture, its toxicity to the cells remains a concern. Hence, using optimum concentrations of components in the freezing medium containing DMSO is imperative and critical.

Abbreviations

AGG	Aggrecan
ALP	Alkaline phosphatase
COLI	Collagen type I
CTGF	Connective tissue growth factor
DEC	Decorin
DM	Diabetes mellitus
ECM	Extracellular matrix
DTA	Ethylenediaminetetraacetic acid
ELA	Elastin
PDGFRa-positive cells	Platelet-derived growth factor receptor-a-positive cells
MSCs	Telomerase transformed mesenchymal stromal cells
AA	L-Ascorbic acid
MSCs	Telomerase transformed mesenchymal stromal cells
TIMP1	Tissue inhibitor of metalloproteinases 1
HAS3	Hyaluronic acid synthase 3

Supplementary Information

The online version contains supplementary material available at https://doi. org/10.1186/s13104-023-06549-y.

Additional file 1: Western Blot Images.

Acknowledgements

The authors would like to thank Mr. Manju N. J. (Sharjah Institute of Medical Research, University of Sharjah) for providing technical support in flow cytometry analysis.

Author contributions

Conceptualization: AS, BMM; methodology: AS, and SHA; software: AS; validation: BMM; writing—original draft preparation: AS, and RR; writing—review and editing: BMM, AAK and SHA; funding acquisition: BMM. All authors have read and agreed to the published version of the manuscript.

Funding

This research was funded by Boehringer Ingelheim (UOS-B12018), University of Sharjah (1701090129-P) and Al Jalila foundation AJF2018082).

Availability of data and materials

All data presented from published data have been subjected to prior approval from the concerned authors.

Declarations

Ethics approval and consent to participate Not applicable.

Consent for publication

Not applicable.

Competing interests

The authors declare no conflict of interest. The funders had no role in the design of the study; in the collection, analyses, or interpretation of data; in the writing of the manuscript, or in the decision to publish the results.

Author details

¹College of Medicine, Research Institute of Medical and Health Sciences, University of Sharjah, P.O. Box: 27272, Sharjah, United Arab Emirates. ²Research Institute of Sciences and Engineering, University of Sharjah, P.O. Box: 27272, Sharjah, United Arab Emirates. ³The Department of Biomedical and Clinical Sciences, Linkoping University, Linkoping 58183, Sweden. ⁴Department of Applied Biology, College of Science, University of Sharjah, P.O. Box: 27272, Sharjah, United Arab Emirates. ⁵BioGrad, 61 Stephenson Way, Liverpool L13 1HN, UK. ⁶Basic Medical Science Department, College of Medicine, University of Sharjah, P.O. Box: 27272, Sharjah, United Arab Emirates.

Received: 3 January 2023 Accepted: 1 October 2023 Published online: 19 October 2023

References

- Mussa BM, Khan AA, Srivastava A, Abdallah SH. Differentiated PDGFRapositive cells: a novel in-vitro model for functional studies of neuronal nitric oxide synthase. Intl J Mol Sci. 2021;22:3514.
- Tong Z, Sant S, Khademhosseini A, Jia X. Controlling the fibroblastic differentiation of mesenchymal stem cells via the combination of fibrous scaffolds and connective tissue growth factor. Tissue Eng Part A. 2011;17:2773–85.
- 3. Camilleri M. Diabetic gastroparesis. N Engl J Med. 2007;356:820-9.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

