

RESEARCH NOTE

Open Access



# Discrimination between some *Mycoplasma* spp. and *Acholeplasma laidlawii* in bovine milk using high resolution melting curve analysis

Abd Al-Bar Al-Farha<sup>1,2\*</sup> , Kiro Petrovski<sup>1,3,4</sup>, Razi Jozani<sup>5</sup>, Andrew Hoare<sup>6</sup> and Farhid Hemmatzadeh<sup>1,3</sup>

## Abstract

**Objectives:** This study aimed to provide a rapid, accurate and cost-effective diagnostic real time polymerase chain reaction-high resolution melting curve assay (PCR-HRM) to identify and distinguish between four different mycoplasmas and *Acholeplasma laidlawii* isolated at cow-level from a single commercial dairy farm in South Australia. One set of genus-level universal primers was designed targeting the 16S ribosomal RNA gene.

**Results:** Real time PCR-HRM analysis was able to identify and distinguish between five different mollicutes, namely *A. laidlawii*, *M. arginini*, *M. bovirhinis*, *M. bovis* and uncultured *Mycoplasma*. Results were confirmed through sequencing. Our developed assay provides rapid and accurate screening for *Mycoplasma* mastitis detection.

**Keywords:** *Mycoplasma*, Mastitis, Cattle, Milk, *Acholeplasma*

## Introduction

*Mycoplasma* mastitis is of emerging significance worldwide, posing significant economic impacts on the dairy industry. Early detection of *Mycoplasma* mastitis is important to disease control strategies [1]. Several *Mycoplasma* spp. are mainly responsible for mastitis, including *M. bovis*, *M. bovoculi*, *M. alkalescence*, *M. alvi*, *M. bovirhinis*, *Mycoplasma species bovine group 7*, *M. californicum*, *M. dispar*, *M. canis*, *M. verecundum*, *M. canadense* and *M. mycoides* subsp. *mycoides* [2]. *Acholeplasma* spp. may be isolated from milk either as a contaminant [3] or as a co-invader with other mycoplasmas [4, 5]. Conventional microbial culture of mollicutes can be laborious and time-consuming with a variety of species-specific growth requirements [6]. Misdiagnosis of *Mycoplasma* using serological detection is common due to the lag period required for antibody

formation. Therefore, a rapid and accurate diagnostic assay is required for screening of mycoplasma in dairy herds. High resolution melting curve analysis (HRM) has been recently developed and widely used for phenotyping at strain or species-level of various organisms including mycoplasmas [7, 8]. However, field isolates of mastitis related mycoplasmas and other milk environmental mollicutes have not been assessed previously using this method.

The aim of this study was to provide a suitable diagnostic real time polymerase chain reaction-high resolution melting curve analysis (PCR-HRM) to identify and distinguish between five different mollicutes isolated at cow-level from a single commercial dairy farms in South Australia.

## Main text

### Methods

Samples were selected based on conventional PCR findings of a previous study conducted on single commercial dairy farm in South Australia. This farm had a history of repeated mastitis treatment failure with high somatic cell

\*Correspondence: abdal-barahmednoori.al-farha@adelaide.edu.au

<sup>1</sup> Present Address: School of Animal and Veterinary Sciences, Roseworthy Campus, The University of Adelaide, Roseworthy, South Australia 5371, Australia

Full list of author information is available at the end of the article

count (SCC) and poor response to antimicrobials [4]. Six isolates for each of the following species were selected in this study: *A. laidlawii*, *M. arginini*, *M. bovirhinis*, *M. bovis* and uncultured mollicutes. One set of genus-level universal primers targeting the 16S rRNA gene was designed for real time-PCR. Forward primer Mol-F: GGC GAAYGGGTGAGTAACAC and reverse primer Mol-R: CATHGYCTTGGTRRGCYNTTA. The real time PCR mixture was prepared using HRM kit AccuMelt HRM SuperMix (Quantabio, Australia). DNA amplification was conducted in a 96 microplate (Illumina, San Diego, CA, USA). Each well contained 10  $\mu$ L reaction solution of 5  $\mu$ L HRM SuperMix, 1  $\mu$ L DNA template (approximately 20 ng), 1  $\mu$ L each primer (0.2 nmol) and 2  $\mu$ L nuclease free water (Qiagen, Germany). The reaction was conducted using an Illumina Thermal Cycler with pre-heating activation for 2 min followed by 40 PCR cycles of three steps: denaturation at 95  $^{\circ}$ C for 15 s, annealing at 60  $^{\circ}$ C for 45 s, then extension at 72  $^{\circ}$ C for 15 s. HRM was performed at 55–95  $^{\circ}$ C at the rate of 0.1  $^{\circ}$ C. Results were analysed via EcoStudy software (version 5.0, Illumina). PCR products were subject to electrophoresis in 1.5% agarose gels and visualised by staining with Gel Red. PCR products from the 16S rRNA gene were submitted to the Australian Genome Research Facility Ltd (AGRF, Adelaide, South Australia) for Sanger sequencing. Each fragment was sequenced in forward and reverse directions. To reconstitute the sequence, forward and reverse sequences were edited and assembled using BioEdit


package v.7.0.4.1. Edited sequences were blasted against existing sequences in GeneBank using the basic local alignment search tool (BLAST) (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>) and nucleotide sequences from relevant *Mycoplasma* strains were used as reference strains for nucleotide alignments using ClustalW program version 2.

## Results

Five different mollicutes, *A. laidlawii*, *M. arginini*, *M. bovirhinis*, *M. bovis* and uncultured mollicutes, produced normalised and derivative melt curves (Fig. 1). *A. laidlawii* (Accession No. LC201977.1) generated one melting peak at 81.2  $^{\circ}$ C, *M. arginini* (Accession No. LC158832.1) generated two melting peaks at 88.5 and 84.7  $^{\circ}$ C. *M. bovirhinis* (Accession No. AP018135.1) generated three melting peaks at 85.7, 77.6 and 88.2  $^{\circ}$ C. *M. bovis* (Accession No. KX462439.1) generated two melting peaks at 77.6 and 85.2  $^{\circ}$ C. Uncultured *Mycoplasma* spp. (Accession No. LT679634.1) generated one melting peak at 83.9  $^{\circ}$ C.

## Discussion

While conventional culture, the traditional method for mollicute detection, imposes technical challenges in distinguishing between milk pathogenic and saprophytic mollicutes, our study indicated that real time PCR-HRM assay provides a sensitive, rapid and cost-effective screening method to identify and discriminate between



some pathogenic and environmental mollicutes isolated from milk DNA. Effects of some of these mollicutes on SCC and milk composition have been previously studied [4]. We considered SCC as the crucial factor that discriminates between contagious and environmental mollicutes. *M. bovis* was widely reported as a main mastitis causing *Mycoplasma* [9, 10]. Inconsistent results have been reported regarding involvement of *M. arginini*, *M. bovirhinis* and *A. laidlawii* in bovine mastitis, particularly with co-infection *Mycoplasma* mastitis [5, 11, 12]. However, several studies indicate these mollicutes are not significant pathogens [3, 13].

Melting profile, introduced in 2002, is widely used for genotyping a wide range of microorganisms [14–17]. HRM-based assay describes correlation between temperature and DNA extent of denaturation [18]. The variety of melting temperatures for different species is attributed to DNA length, sequencing and GC content [14]. In summary, as an alternative to sequencing, our developed real time PCR-HRM assay offered a rapid, low-cost and simple discriminative method to distinguish between some mastitis causing pathogenic mycoplasmas and other saprophytic mollicutes in bovine milk. This method was useful for screening of *Mycoplasma* mastitis and can be extended to identify more mollicutes species.

### Limitations

One of the limitations of HRM-based analysis in *Mycoplasma* mastitis detection is the inability to detect co-infection cases due to amplicon concentration differences and the requirement of separation each individual amplicon Tm [18]. Primers used in this study were designed to target more spp. of major *Mycoplasma*-causing mastitis in dairy herds. However, in this study, we used only field isolates of five different mollicutes.

### Abbreviations

HRM: high resolution melt; PCR: polymerase chain reaction; SCC: somatic cell count.

### Authors' contributions

AAA, FH and KP participated in the study design and coordination. AAA, KP and AH contributed to sample collection. AAA and FH contributed in sample processing and real time PCR in the PC2 laboratory. RJ designed the primers for HRM analysis. All authors were involved in drafting the manuscript, corrected. All authors read and approved the final manuscript.

### Author details

<sup>1</sup> Present Address: School of Animal and Veterinary Sciences, Roseworthy Campus, The University of Adelaide, Roseworthy, South Australia 5371, Australia. <sup>2</sup> Mosul Technical Institute, Technical Foundation, Northern Technical University, Mosul 41000, Iraq. <sup>3</sup> Australian Centre for Antimicrobial Resistance Ecology, The University of Adelaide, Adelaide, South Australia 5000, Australia. <sup>4</sup> Davies Centre, School of Animal and Veterinary Sciences, The University of Adelaide, Roseworthy, South Australia 5371, Australia. <sup>5</sup> Department of Veterinary Clinical Sciences, University of Tabriz, Tabriz 51666-14766, Iran. <sup>6</sup> South East Vets, Mt Gambier, South Australia 5290, Australia.

### Acknowledgements

The authors would like to express their special thanks to the farmer and farm staff for their cooperation and interest in the study.

### Competing interests

The authors declare that they have no competing interests.

### Availability of data and materials

The datasets generated during and/or analysed during the current study are available from the corresponding author on reasonable request.

### Consent to publish

Not applicable.

### Ethics approval and consent to participate

Not applicable. Samples were collected by field veterinary support as a part of the mastitis investigation as per farmer request and no Animal Ethics application was required (Australian code for the care and use of animals for scientific purposes, 8th edition, 2013).

### Funding

The first author was partially funded by The Higher Committee for Education Development in Iraq in form of a scholarship. All field and laboratory procedures were financially supported by the University of Adelaide.

### Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 30 December 2017 Accepted: 1 February 2018

Published online: 07 February 2018

### References

- Kirk JH, Lauerman LH, Roberts C. *Mycoplasma* mastitis in dairy cows. *Compend Contin Educ Pract Vet.* 1994;16(4):541–52.
- Nicholas R, Ayling R, McAuliffe L. *Mycoplasma* mastitis. *Vet Rec.* 2007;160(11):382 (author reply 383).
- Boonyayatra S, Fox LK, Gay JM, Sawant A, Besser TE. Discrimination between *Mycoplasma* and *Acholeplasma* species of bovine origin using digitonin disc diffusion assay, nisin disc diffusion assay, and conventional polymerase chain reaction. *J Vet Diagn Investig.* 2012;24(1):7–13.
- Al-Farha AAB, Hemmatzadeh F, Khazandi M, Hoare A, Petrovski K. Evaluation of effects of *Mycoplasma* mastitis on milk composition in dairy cattle from South Australia. *BMC Vet Res.* 2017;13(1):351.
- Counter DE. A severe outbreak of bovine mastitis associated with *Mycoplasma bovis* and *Acholeplasma laidlawii*. *Vet Rec.* 1978;103(7):130–1.
- Hahn RG, Kenny GE. Differences in arginine requirement for growth among arginine-utilizing *Mycoplasma* species. *J Bacteriol.* 1974;117(2):611–8.
- Ghorashi SA, Noormohammadi AH, Markham PF. Differentiation of *Mycoplasma gallisepticum* strains using PCR and high-resolution melting curve analysis. *Microbiology.* 2010;156(Pt 4):1019–29.
- Rebelo AR, Parker L, Cai HY. Use of high-resolution melting curve analysis to identify *Mycoplasma* species commonly isolated from ruminant, avian, and canine samples. *J Vet Diagn Investig.* 2011;23(5):932–6.
- Aebi M, van den Borne BH, Raemy A, Steiner A, Pilo P, Bodmer M. *Mycoplasma bovis* infections in Swiss dairy cattle: a clinical investigation. *Acta Vet Scand.* 2015;57:10.
- Kunkel JR. Isolation of *Mycoplasma bovis* from bulk milk. *Cornell Vet.* 1985;75(3):398–400.
- Szacawa E, Niemczuk K, Dudek K, Bednarek D, Rosales R, Ayling R. *Mycoplasma bovis* infections and co-infections with other *Mycoplasma* spp. with different clinical manifestations in affected cattle herds in eastern region of Poland. *Bull Vet Inst Pulawy.* 2015;59(3):331–8.
- González RN, Wilson DJ. Mycoplasmal mastitis in dairy herds. *Vet Clin N Am Food Anim Pract.* 2003;19(1):199–221.
- Jasper D, Dellinger J, Rollins M, Hakanson H. Prevalence of mycoplasmal bovine mastitis in California. *Am J Vet Res.* 1979;40(7):1043–7.

14. Reed GH, Kent JO, Wittwer CT. High-resolution DNA melting analysis for simple and efficient molecular diagnostics. *Pharmacogenomics*. 2007;8(6):597.
15. Ren X, Fu Y, Xu C, Feng Z, Li M, Zhang L, Zhang J, Liao M. High resolution melting (HRM) analysis as a new tool for rapid identification of *Salmonella enterica* serovar *Gallinarum* biovars *Pullorum* and *Gallinarum*. *Poult Sci*. 2017;96(5):1088–93.
16. Yong TB, Hashim R, Noor AM, Hamzah SH, Ahmad N. Identification of *Brucella* spp. isolated from human brucellosis in Malaysia using high-resolution melt (HRM) analysis. *Diagn Microbiol Infect Dis*. 2015;81(4):227–33.
17. Gago S, Zaragoza Ó, Cuesta I, Rodríguez-Tudela JL, Cuenca-Estrella M, Buitrago MJ. High-resolution melting analysis for identification of the *Cryptococcus neoformans*–*Cryptococcus gattii* complex. *J Clin Microbiol*. 2011;49(10):3663–6.
18. Tong SY, Giffard PM. Microbiological applications of high-resolution melting analysis. *J Clin Microbiol*. 2012;50(11):3418–21.

Submit your next manuscript to BioMed Central  
and we will help you at every step:

- We accept pre-submission inquiries
- Our selector tool helps you to find the most relevant journal
- We provide round the clock customer support
- Convenient online submission
- Thorough peer review
- Inclusion in PubMed and all major indexing services
- Maximum visibility for your research

Submit your manuscript at  
[www.biomedcentral.com/submit](http://www.biomedcentral.com/submit)

